## **User Instruction Manual**

Oxford® Up

To avoid injury, read user manual prior to use.

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
## Contents

1. Introduction .....	3
2. Assembly & Commissioning Instructions.....	4
3. Safety Precautions.....	5
4. Operating Instructions.....	7
5. Maintenance Schedule .....	15
6. Technical Specifications .....	16
7. Warranty .....	17

# 1. Introduction

The Oxford Up is an assistive standing and transfer device, designed for moving a person over short distances. This might be from bed to chair, chair to wheelchair or to and from a toilet/commode, for example. The Up has a compact footprint, enabling effective use, even in confined spaces.

The packing consists of a strong, purpose built carton that is used for both export and domestic markets to ensure the safe arrival of the product.

To properly maintain your product please refer to the maintenance schedule included in this document. If you are at all unsure what your country's servicing requirements are, please check with your authorised distributor and/or a local government agency.

**The CE mark:** The Oxford Up carries the CE mark and complies with the following EC directives:

- Medical Device Directive MDR 2017/745

## Statement of Intended Use

The Oxford Up is intended for use by an individual who has the ability to stand unassisted or with minimal assistance. It will enable their transfer from one resting surface to another, such as from a bed to a chair. It is an active product for users who can contribute effort to the standing and transfer process. By encouraging user participation, the Up is a valuable rehabilitation aid, improving natural function and tone. The Up is designed for indoor use and for transfers on flat surfaces, over short distances only.

## Expected Service Life


Oxford Professional products are designed and tested for a minimum service life of ten (10) years, subject to the use and maintenance procedures stated in this manual. Use, other than in accordance with these instructions, may compromise service life.

## 2. Assembly & Commissioning Instructions

### Carton Contents

Place the carton in a clear working area and open carefully. The carton contains:

- Oxford Up
- User Guide (This document)

### Preparation for Use

#### **⚠ WARNING**

**The Oxford Up will need to be lifted with care. You may need assistance to lift the Oxford Up from its carton.**

The Up requires no additional assembly once removed from its packaging; it is ready for use in accordance with the guidelines and daily checks contained in this document.

### 3. Safety Precautions

Please read and follow the safety precautions listed below. The operation and use of Oxford products is simple and straightforward. Following these few basic safety precautions will make transfer operations easy and trouble free.

**READ AND UNDERSTAND THE USER INSTRUCTIONS BEFORE USING YOUR UP.**

**WARNING:** Important safety information for hazards that might cause serious injury.

**CAUTION:** Information for preventing damage to the product.

**NOTE:** Information to which you should pay special attention.

#### **WARNING**

- **ALWAYS** carry out a full and thorough risk assessment prior to using the Up.
- In order to use the Up safely, the patient must display the following abilities:
  - Support their own body weight while standing
  - Maintain balance for an extended period while standing (with support)
  - A firm grip in order to grasp the frame handles for support and stability
  - Ability to understand and carry out instruction
- The Up is intended for short transfers on flat/even floors.
- **ALWAYS** plan your transfer operations before commencing.
- **ALWAYS** carry out the **DAILY CHECK LIST** before using this product.
- **ALWAYS** familiarise yourself with the operating devices and safety features of this product before use.
- **DO NOT** use a sling or belt unless it is recommended for use with the device.
- **ALWAYS** check the sling/belt is suitable for the particular patient and is of the correct size and capacity.
- **NEVER** use a sling/belt which is frayed or damaged.
- **ALWAYS** fit a sling/belt according to the user instructions provided.
- **ALWAYS** check the safe working load of the product(s) is/are suitable for the weight of the patient.
- **ALWAYS** carry out transfers according to the instructions in the user manual.
- **NEVER** attempt to bypass a safety feature because it seems easier to operate the product.
- **NEVER** force an operating or safety device. All devices are easy to use and do not require excessive force to operate. If a device is not working easily there will be a reason. Forcing will only strain or damage the product and may compromise safety.
- **ALWAYS** apply both brakes after positioning the Up into position in front of the patient prior to asking them to stand.
- When transferring to or from a wheeled/movable surface such as a commode, wheel chair or bed, ensure the brakes (of that surface) are applied to avoid any inadvertent movement and ensure a safer, more stable transfer.


- **DO NOT** attempt to manoeuvre the device by pushing on the patient.
- **ALWAYS** manoeuvre the device using the handle provided. **DO NOT** attempt to apply side force to the base or frame of the device in order to manoeuvre.
- **DO NOT** push a loaded device at speeds exceeding a slow walking pace (3 kilometres/hour 0.8 metres/second).
- **DO NOT** push the device over wet surfaces, uneven or rough ground or surfaces littered with debris, particularly when loaded.
- **DO NOT** attempt to push or pull a loaded device over a floor obstruction of which the castors are unable to ride over easily.
- **DO NOT** bump the device down steps, loaded or unloaded.
- **DO NOT** attempt to negotiate a loaded device on a slope which exceeds 1:12 (approximately 5 degrees). Joerns Healthcare recommend a second helper is present when moving a patient on a slope.
- **DO NOT** park a loaded device on ANY sloping surface.
- **DO NOT** transfer a patient unless you are trained and competent to do so.
- Your device is for patient standing and transfers. **DO NOT** use it, or allow it to be used, for any other purpose.
- **DO NOT** use this product for self transfers; a patient must always be assisted by caregiver(s).
- **ENSURE** the caregiver(s) wear appropriate footwear when using this product.
- **OXFORD RECOMMENDS THE USE OF GENUINE OXFORD PARTS.** Oxford sling/belt and transfer products are designed to be compatible with one another. For country specific guidance on sling/belt use and compatibility, please refer to the sling/belt label or contact your local market distributor or Joerns Healthcare.

## 4. Operating Instructions

### 1. Swing-away Seat Pads


The seat pads are designed to swing-away allowing easy access for the patient when standing. There is an integrated handle on the back of each seat pad to facilitate this.

Once the patient is in position and standing, the seat pads can be swung back into the seated position, allowing the patient to sit down.


### 2. Opening & Closing the Legs

The legs of the Up can be opened and closed to allow easier positioning around furniture. To open the legs, push down the left hand pedal. To close, push down the right hand pedal.

### 3. Applying the Brakes

Each castor can be braked independently by pressing down firmly on the integrated brake pedal. To release the brakes, simply manoeuvre the pedal back up.

#### 4. Approaching the Patient with the Oxford Up

The following guidance applies no matter what the surface you are transferring from, whether it be a chair, wheelchair, bed, toilet or commode.


### **⚠ WARNING**

**ALWAYS** carry out a full and thorough risk assessment prior to using the Up.


Maintain a dialogue with the patient, ensuring they understand what is happening and they are comfortable and at ease.

If a sling has been deemed necessary as part of the risk assessment, ensure you have the appropriate, compatible style and size.


Rotate the swing-away seat pads into the open position to allow patient access.


To gain improved access around large chairs and bulky commodes, open the legs of the Up before approaching the patient.

Carefully wheel the Up into position, in front of the patient, placing their feet on the footplate with their knees resting comfortably against the cushioned knee pad.

Once the patient is comfortably in position, apply both castor brakes.

## 5. Assisted Standing and Transfers with the Oxford Up


The following guidance applies regardless of the surface you are transferring to, whether it be a chair, wheelchair, bed, toilet or commode.

### **⚠ WARNING**


**When transferring to or from a wheeled/movable surface such as a commode, wheelchair or bed, ensure the brakes (of that surface) are applied to avoid any inadvertent movement and ensure a safer, more stable transfer.**

Ask the patient to grasp the crossbar and prepare to raise themselves to a standing position. Offer words of encouragement and provide assistance as required.


**NOTE:** The patient's knees should remain in contact with the knee pad throughout the standing and transfer process.

Once the patient is standing, rotate the swing-away seat pads into the seated position before asking the patient to carefully lower themselves into a seated position.


Ensure the patient is comfortable and fully supported before releasing the brakes and proceeding carefully with the transfer.


### **⚠ WARNING**

**When transferring to or from a wheeled/movable surface such as a commode, wheelchair or bed, ensure the brakes (of that surface) are applied to avoid any inadvertent movement and ensure a safer, more stable transfer.**


Wheel the patient directly in front of the target surface and apply the brakes. Open the Up's legs to help navigate and get close to wider furniture including bulky commodes and armchairs.

Ask the patient to grasp the crossbar and raise themselves into a standing position, providing assistance if required.


Once the patient is standing, rotate the swing-away seat pads into the up position, providing a barrier-free space in behind.

**⚠ WARNING**


**When transferring to or from a wheeled/movable surface such as a commode, wheelchair or bed, ensure the brakes (of that surface) are applied to avoid any inadvertent movement and ensure a safer, more stable transfer.**

With the patient continuing to grasp the crossbar, ask them to carefully lower themselves on to the target surface, providing assistance as necessary.

## 6. Sling Compatibility & Attachment

The **Oxford Deluxe Standing sling with clips** is suitable for use with the Oxford Up. Various sizes are available depending on patient requirements.

### ⚠ WARNING

Joerns Healthcare recommends that slings be checked regularly and particularly before use for signs of fraying or damage. **DO NOT** use slings that are worn or damaged.

### ⚠ WARNING

**OXFORD RECOMMENDS THE USE OF GENUINE OXFORD PARTS.** Oxford sling, lift and transfer products are designed to be compatible with one another. For country specific guidance on sling use and compatibility, please refer to the sling label or contact your local market distributor or Joerns Healthcare.

### ⚠ WARNING


Refer to maximum weight capacity of product. Sling capacity is limited by the maximum capacity of the product it is used on.

### Fitting Instructions:


**NOTE:** For full fitting instructions, please refer to the user guide supplied with each sling.

If your patient would benefit from additional reassurance and support, the Oxford Deluxe Standing sling with clips can be used in combination with the Up.

Place the sling around the patient just above the waist line and secure the waist belt around the patient using the clip. This strap need not be tight. It is there to hold the sling in position whilst attaching the sling to the Up. Ensure that the client's arms are outside the sling.

When fitted to the patient, the sling is connected to the clip attachment pins by using the Securi3 clip system (step by step guide as follows). To remove, simply reverse this procedure.


SECURI3 CLIP SYSTEM

**PLEASE NOTE:**

1. Attach Securi3 clips to the attachment pins as described in the illustration shown.
2. Ensure all green safety loops are visible before transferring.

**NOTE:** Use the appropriate clip to achieve a natural seated position that provides good, comfortable support without pulling the patient forward or allowing them to lean back and become unbalanced.

Continue with the transfer according to the guidance in this document.

## 5. Maintenance Schedule & Daily Check List

All Oxford products are designed for minimum maintenance, however some safety checks and procedures are required. A schedule of DAILY tasks are detailed below.

Daily checks and regular inspection will ensure your product is kept in optimum safe working condition. A list of spare parts is available upon request.

### Daily Check List:

Joerns Healthcare strongly recommends the following checks be carried out on a daily basis and before using the Up.

- **ENSURE** the product moves freely on its castors.
- **ENSURE** both brakes are operational and engage correctly.
- **ENSURE** the legs open and close correctly.
- **EXAMINE** the sling attachment points for excessive wear. If in doubt - do not use.
- **EXAMINE** slings for fraying or other damage. DO NOT use any sling if damaged or if the sling shows signs of wear.
- **INSPECT** the frame and base for signs of cracks, damage or distortion.

### Spare Parts & Repairs:

#### **WARNING**

**JOERNS HEALTHCARE RECOMMENDS THE USE OF GENUINE OXFORD PARTS.**

Please contact your local authorised Oxford service provider to arrange replacement parts and repairs. If you are unsure of the organisation you need to contact, our customer service team will be more than happy to assist you. Contact details can be found on the reverse cover of this manual.

### Cleaning:

Using a clean, damp cloth, wipe clean with ordinary soap and water and/or any hard surface disinfectant. Harsh chemical cleaners or abrasives should be avoided as these may damage the surface finish of the product.

## 6. Technical Specifications

Safe Working Load.....	200kgs
Overall Width.....	640mm
Overall Height .....	1065mm
Overall Length.....	955mm
Legs Open - External Width.....	850mm
Legs Open - Internal Width.....	730mm
Legs Closed - External Width.....	480mm
Legs Closed - Internal Width.....	350mm
Leg Height.....	95mm
Ground Clearance.....	15mm
Overall Height Of Base (Footplate).....	65mm
Knee Pad Height (To top of knee pad).....	465mm
Turning Radius.....	1165mm
Front Castors.....	75mm
Rear Castors .....	100mm

### Weights:

Total.....	22.5kgs
------------	---------

### Environmental Conditions:

Outside this environment functionality and safety may be compromised.

### Operating:

Temperature .....	5°C to 40°C
Relative humidity.....	20% to 90% @ 30°C - not condensing

### Storage:

Temperature .....	-20°C to +60°C
Relative humidity.....	20% to 90% @ 30°C - not condensing

## 7. Warranty

Joerns Healthcare has an established network of authorised distributors and service agents who will be pleased to handle all of your purchasing, warranty, repair and maintenance enquiries. It is recommended that our products are commissioned by your authorised distributor or service agent.

The distributor or service agent operates the warranty programme, so it is important to keep a record of their name, address and telephone number so they can be contacted should any problem arise. If you are in any doubt as to where your product was purchased, Joerns Healthcare can trace the supplier if you quote the serial number of the product.

All warranty claims are subject to the correct use and maintenance in accordance with the user instruction manual provided with each product. Damage caused by use in unsuitable environmental conditions or failure to maintain the product in accordance with user and service instructions is not covered. Your statutory rights remain unaffected.

**REMEMBER:** Contact your authorised distributor or service agent for purchases, repairs, servicing and certified maintenance.

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