

User Manual - Manuel d'utilisation  
Gebrauchsanweisung - Gebruikershandleiding  
Användarhandbok - Brugermanual  
Manuale d'uso - Manual de Usuario  
Manual do Utilizador


**Primo**  
**P02033 / P02034 / P02044 / P02047**

**170924**  
Rev.002

# User Manual

## Primo

P02033 /P02034 /P02044 /P02047

Hill-Rom Industries S.A.  
188, Rue du Caducée  
Parc Euromédecine  
34195 MONTPELLIER Cedex 5 - FRANCE  
Tel: +33 (0)4 67 04 64 04  
Fax: +33 (0)4 67 04 64 00  
[www.hill-rom.com](http://www.hill-rom.com)

First edition, first printing: june 2012  
170924, Rev.002 / October 2012

The information contained in this manual is confidential and may not be reproduced or divulged in any form or by any means without the prior written permission of Hill-Rom.

Hill-Rom® is a registered trademark of Hill-Rom Services, Inc.

I-mmersion™ is a trademark of Hill-Rom Services, Inc.

Primo™ is a trademark of Hill-Rom Services, Inc.

Velcro® is a registered trademark of Velcro Industries BV.

Hill-Rom reserves the right to make changes to the design, characteristics and models without prior notice. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

To order copies of this manual, refer to the last page, identify your national Hill-Rom representative and order the article with the part number 170924.

© 2012 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

# Table of Contents

## Introduction, specifications

Structure of the manual . . . . .	1
Symbol definitions . . . . .	2
Applications . . . . .	3
Description . . . . .	3
Indication . . . . .	3
Contraindication . . . . .	3
Safety and Usage Tips . . . . .	4
Putting into service . . . . .	4
Risk prevention . . . . .	4
Electrical safety . . . . .	5
Technical specifications . . . . .	6
Essential performances . . . . .	6
Therapeutic Mattress . . . . .	6
Control unit . . . . .	7
Conditions for transport, storage and use . . . . .	8
Special conditions of use . . . . .	8
Product overview . . . . .	9
General symbols on the device . . . . .	10
General symbols on the upper cover . . . . .	10
General symbols on the control unit . . . . .	10
Buttons and light indicators on the control unit . . . . .	11
General symbols on the labels . . . . .	12
Product identification labels . . . . .	13
Control unit . . . . .	13
Packaging . . . . .	13
Upper and bottom covers . . . . .	14

## Installing the patient

Before placing the patient on the Primo™ . . . . .	15
List of compatible bed frame . . . . .	15
Accessories . . . . .	15
Installation . . . . .	16

## Mobilizing and securing the patient

Functions . . . . .	19
Therapeutic mode . . . . .	19
Maximum inflation . . . . .	19
CPR . . . . .	20
Alarms . . . . .	21
Alarm silence . . . . .	21
Mains power fault . . . . .	21
Malfunction . . . . .	21


## Moving the device

Transport mode . . . . .	23
--------------------------	----

Transferring the mattress . . . . .	23
Storing the mattress . . . . .	24
<b>Disinfection, servicing</b>	
Cleaning and Disinfection . . . . .	25
Safety Tips . . . . .	25
Infection control . . . . .	25
General recommendations . . . . .	26
Frequency of cleaning and disinfection of different parts of the mattress . . . . .	27
Clean and disinfect after the departure of a patient or when changing the linen. . . . .	27
Cleaning and disinfection at regular intervals or in the event of high risk of contamination. . . . .	28
Final steps . . . . .	30
Maintenance . . . . .	30
Safety recommendations . . . . .	30
Preventive maintenance. . . . .	30
Troubleshooting . . . . .	31
De-commissioning . . . . .	32
<b>Accessories and spare parts</b>	
Seat cushion . . . . .	33
Installation . . . . .	33
Deflating and packing the seat cushion . . . . .	34
Cleaning and Disinfection. . . . .	34
Upper covers. . . . .	35
Quick inflate/deflate kit . . . . .	35
Transport bag . . . . .	36
<b>Appendix</b>	
Warranty and after sales service conditions . . . . .	37
Regulatory requirements . . . . .	38
Manufacturer's declaration - Electromagnetic Compatibility . . . . .	39
Electromagnetic emissions compliance . . . . .	39
Electromagnetic immunity compliance. . . . .	39
Recommended separation distances . . . . .	40

## Structure of the manual


For every type of use, Hill-Rom mattresses provide patients with optimal comfort and greater independence for a feeling of well-being that is conducive to a swift recovery. They are also easy to use for caregivers.

## Symbol definitions

This manual contains different typefaces and icons designed to improve readability and increase understanding of its content. Note the following examples:

- Standard text: used for regular information.
- **Bold text:** emphasizes a word or phrase.
- ⓘ highlights special information or explains very important instructions.

The symbols below represent different risks or hazards:

Symbol	Description
	<b>WARNING</b> <ul style="list-style-type: none"><li>• This symbol indicates that the failure to follow the associated recommendation can put the patient or the user in danger, or damage the equipment.</li></ul>
	<b>CAUTION</b> <ul style="list-style-type: none"><li>• This symbol indicates that the failure to follow the associated recommendation can result in damage to the equipment.</li></ul>
	Tip
	Risk of falling
	Caught hazard warning
	Risk of crushing an upper limb
	Chemical Hazard Warning
	Electric Shock Hazard
	Biological danger

---

# Applications

---

## Description

The Primo™ is a therapeutic mattress equipped with a pressure control system designed to operate in continuous (or static) low-pressure mode. The Primo™ is a therapeutic solution by Immersion™ that distributes pressure and reduces the shearing effects on the surface of the patient's skin.

The Primo™ comprises a mattress and an external control unit.

It offers the following functionality:

- a CPR valve to facilitate external cardiac massages;
- a maximum inflation mode for greater patient mobility.

---

## Indication

The Primo™ therapeutic mattress helps to prevent and treat pressure ulcers in low to moderate risk adult patients, within the recommended patient weight limits of 30 to 150 kg, in order to achieve clinical performance assessed in all the usual positions of the adjustable head section.

It can be used as a mattress in the following environments, as defined in the standard IEC 60601-2-52:

- application environment 2 (short-term care in hospitals or other medical establishments);
- application environment 3 (long-term care in medical establishments);
- application environment 4 (care administered at home);
- application environment 5 (outpatients or ambulatory care).

The Primo™ is intended to be used with a sheet between the patient skin and the surface of the mattress.

**i** *In accordance with the NPUAP/EPUAP directives<sup>1</sup>, we recommend that the condition of each patient be regularly checked. For patients with special needs, we recommend the use of the more suitable Immersion™ Therapy system. Caregivers are responsible for taking this decision, in accordance with modern care practices.*

---

## Contraindication

**The Primo™ must not be used with patients:**

- suffering from an unstable fracture of the spine.  
For any other unstable fractures, a medical examination is necessary to decide whether the use of the device is appropriate;
- with atypical anatomies;
- suffering from cervical or trans-bone traction.

---

<sup>1</sup> NPUAP / EPUAP - Pressure Ulcer Prevention - Quick Reference Guide, January 2010  
NPUAP / EPUAP - Pressure Ulcer Treatment - Quick Reference Guide, December 2009

---

# Safety and Usage Tips

---

## Putting into service

Before using the mattress, it is essential to have a thorough understanding of this manual. This manual contains instructions for general use and maintenance and guarantees improved safety. Caregivers must have access to this manual.

Caregivers must be informed of the risks that may be encountered in the use of electric devices. Product training can be provided on demand.


When using the mattress with medical devices (accessories), the user must ensure that safety and conformity requirements are met.

Before using the mattress for the first time or when removing from storage:

- check the condition and conformity of the electrical system with the applicable safety standards;
- connect the control unit of the Primo™ device to the mains (See “Electrical safety” page 5).
- make sure that all the functions of the device (See “Installation” page 16) are in good conditions of operation;
- make sure that the device and the care environment are in a good state of hygiene (See “Disinfection” page 28).

---

## Risk prevention

***Improper use of the device Primo™ can result in risks for the patient or the user. The following recommendations must be read and followed.***

In view of the multitude of models of frames and siderails, and for safety reasons, Hill-Rom advises that all necessary precautions must be taken, especially with regard to the height of the siderails and the dimensions of the sleep surface. If this device is used on a bed with siderails that are less than 22cm above the mattress, patients must not be left unattended.

For safety reasons, it is advisable to use the lock-out functions of the bed in the following situations:

- during all interventions on the patient or the device (e.g.: examinations, transfers, maintenance);
- if the patient suffers from particular behavioral difficulties (e.g. : agitation, confusion, disorientation, obsessional behavior and aged and frail patients).

Properly trained medical staff should determine how the device should be used and the required level of monitoring or constraint.

It is imperative to observe the practices relating to the safety of caregivers. Special attention must be taken when redistributing the load application points, since there is a danger that the bed will tip when the frame is moved.

Surface impermeability and therapy capabilities could be affected by needle sticks or other mattress bladder punctures. Caregivers should be instructed to prevent mattress bladder punctures caused by needle sticks.


*Avoid mechanical impacts.*

## Electrical safety


The mains power supply must comply with following standards:

- NF C15-100 and NF C15-211 (France);
- International Electrotechnical Commission (IEC) 60364 for other countries.

Check that the power requirements of the device Primo™ on the identification label correspond to the power supply voltage of the hospital (See “Product identification labels” page 13).

The Primo™ should be connected to a power system equipped with a maximum 30 mA earth leakage circuit breaker, in compliance with IEC 60364-5-53.

*Connect the device to the nearest wall plug in order to leave the shortest possible length of cable on the floor.*

In compliance with standards relating to electromagnetic interference for medical devices, this product does not interfere with other medical devices or is not susceptible to interference when combined with other medical devices that also comply with the electromagnetic standards in force.

Some devices, particularly older ones that do not comply with the electromagnetic compatibility standards, may however undergo interference or may themselves interfere with the working of this product.

The users of such devices are responsible for ensuring that any malfunctions will not endanger the patient or any other person.

When intravascular or intracardiac connections are in use, the electric potentials of all the unprotected metal parts of the appliance and the bed need to be equalized.

This label indicates that oxygen tents must never be used and that only the use of nasal tubes and oxygen masks is authorized. For reasons of safety, masks and tubes should always be kept at a higher level than the mattress support platform.

## Technical specifications

**i** Hill-Rom has a continuous improvement policy for its products. Therefore specifications are liable to be altered without notice.


## Essential performances

Percentile	P5-P30* Woman		P75-P95* Man	
Angle of the adjustable head section	0°	45°	0°	45°
I-mmersion™ setting	22 mbar	26 mbar	35 mbar	40 mbar

\* Data source : IFTH 2006 National Body Measurement Campaign in France (target population: men and women aged between 18 and 70).


**i** Values in mbar +/- 8%.

## Therapeutic Mattress

Characteristics	Description			
	P02033	P02047	P02034	P02044
Length (inflated)	200 cm / 79"			200 cm / 79"
Width (inflated)	85 cm / 33.5"			77 cm / 30.5"
Height (inflated)	16 cm / 6"			16 cm / 6"
Weight	8.4 kg / 18.5 lb			
Autonomy in transport mode	2 hours			
Upper cover	Polyurethane coating on polyester material Low-friction, stretchable in both directions, breathing, bacteriostatic, fungistatic and antimicrobial. Can be wiped and washed.			
Bladders	Polyurethane			
 = 250Kg	Safe working load, including the total weight of the patient, accessories (if they are supported by the support system of the medical bed) and the load supported by these accessories (excluding the weight of the patient).			
Maximum pressure of relief valve	1psi / 69 mbar			
Degree of protection against electric shock	Type BF applied parts protected against defibrillation shocks (labeled 1 and 6 on page 9).			
Degree of protection provided by the cover (IEC 60529)	IP24: protected against access to hazardous parts with a finger and against splashing water.			

## Control unit


### Front view of the control unit of the Primo™

Characteristics	Description			
Model	P02033	P02047	P02034	P02044
Dimensions	15x29x12 cm / 6x11.5x5"			
Weight	3.5 kg / 8 lb			
Power supply	220-240V 50/60 HZ	220-240V 50/60HZ	120V 60HZ	220-240V 50/60HZ
Apparent power	20 VA	20 VA	32 VA	20 VA
Maximum energy consumed	16.4 Wh	16.4 Wh	16.4 Wh	16.4 Wh
Operation of the device	Continuous			
Sleeve material	ABS PC V0			
Device : sound pressure level / sound power (ISO 3744)	41 dB(A) / 52 dB(A)			
Alarm : sound pressure level / sound power (IEC60601-1-8)	56 dB(A) / 67 dB(A)			
Fuses	T160 mA	T160 mA	T 200mA	T160 mA
Compression	0 - 180 mbar	0 - 180 mbar	0 - 215 mbar	0 - 180 mbar
Compressor flow rate	10 l/min	10 l/min	12 l/min	10 l/min
Device classification according to IEC 60601-1	Class II			
Degree of protection provided by the cover (IEC 60529)	IP21 : protected against access to hazardous parts with a finger and against vertically falling water drops.			
Protection against inflammable anesthetic mixtures	Not for use with flammable anesthetics.			
Battery life	1 hour			

---

## Conditions for transport, storage and use

Symbol	Features	Use	Transport, storage <sup>a</sup>
	Temperature	+5°C - +40°C	-25°C - +70°C
	Humidity	15% - 93%	0% - 93 %
	Atmospheric pressure	700 mbar - 1,060 mbar	700 mbar - 1,060 mbar

a. Applicable only if the device is stored in its original packaging.

The device is designed for indoor use only.

*The device must be stored in its original packaging:*


- *protected against light and damp,*
- *at least 10 cm above floor level to prevent fluid ingress;*
- *protected against dust;*
- *outside passageways.*

*Never stack more than 5 mattresses.*

---

## Special conditions of use

A number of special precautions must be taken when using the device in a domestic environment.

*Strangulation warning:*

*Never leave children aged under 3 or animals alone in the vicinity of the device. They may become caught up in the power cord and be strangled.*

*Suffocation warning:*

*Never leave the electric connection accessories supplied with the device within the reach of children or animals.*


*Do not use the device in dusty environments and avoid the presence of soft toys and ash (cigarettes, etc.).*

*Do not expose the device to excessive light, including sunlight.*

*Do not install the device close to a source of intense heat (fireplace, heaters, etc.) or of excessive humidity (sprayers, kettles, etc.).*

*Do not install the device close to pets, insects or pests.  
Protect the device against misuse by children.*



## Product overview


Item	Name	Item	Name
1	Upper cover	5	Foam sub-mattress
2	Air mattress	6	Bottom cover
3	Control unit	7	CPR valve
4	Sensor	8	Straps
		9	Sensor box


## General symbols on the device

### General symbols on the upper cover

	CPR symbol
	Do not walk on or run over the power cord
	Adjust the head and foot end straps
	Foot end
	Zone for notes
	Always install the mattress so that the Immersion™ Therapy text is visible


### General symbols on the control unit

A		Connection for inflation of the seat cushion
B		Fuse 220-240 VAC
		Fuse 120 VAC

## Buttons and light indicators on the control unit

The control unit is very easy to use and, thanks to the automatic pressure regulation, no manual adjustments are necessary

### Buttons

1 - Alarm silence button :  
Inhibits the audible alarms (see “Alarm silence” page 21).

2 - On/Off Button .

### NOTE:

*To switch off the control unit:*

*Press  to silence the battery, in this case no mains power fault alarm will be reported ;*

*Unplug the control unit from the mains.*

3 - Maximum inflation button :  
Increases the pressure inside the mattress to the maximum value for 5 minutes. See “Maximum inflation” page 19.  
Is also used to inflate the seat cushion (see “Seat cushion” page 33).

### Indicator lights / alarms

A - Malfunction alarm :  
In case of malfunction, the audible alarm sounds and the visual alarm yellow comes on (see “Malfunction” page 21).

B - Mains power fault alarm :  
In case of power failure, the audible alarm sounds and the visual alarm yellow comes on (see “Mains power fault” page 21).


C - On/Off Indicator:

The green indicator light at the bottom comes on when the control unit is on.

D - Inflation Indicator :

The green indicator light flashes for 5 minutes while the inflation.

## General symbols on the labels

	Manufacturer		Control unit protected against access to hazardous parts with a finger and against vertically falling water drops.
	Product reference <sup>a</sup>		Cover protected against access to hazardous parts with a finger and against splashing water.
	Serial number		Type BF applied part protected against defibrillation shocks (mattress).
	Alternating Current		Class II Device
	CAUTION, read the safety instructions carefully		Compliant with Directive 93/42/EEC
	Refer to the user manual.		Temperature range
	DO NOT DISCARD. Obey local recycling rules.		Atmospheric pressure range
	Interior use		Humidity range
	Warning : high voltage ! (label on the power supply card in the control unit)		Width of the mattress (model P02044 only)


- a. The product part number contains the following information:
- P+5 figures = model,
  - B = product version letter.

 See "Product identification labels" page 13.

## Product identification labels


### Control unit

To identify your device model and serial number, refer to the identification label on the rear of the control unit.

The label showing the conditions of use and specifications of the device is also located on the rear of the control unit (right-hand side).

See “General symbols on the labels” page 12.

### Packaging

To check the model of the device on the packaging label, refer to the example below:

See symbols in paragraph “Conditions for transport, storage and use” page 8.

## Upper and bottom covers

Open the zip fastener on the mattress to access the upper and bottom covers.

Upper Cover		
P02033/P02034	P02044	P02047
<p><b>Hill-Rom Industries S.A</b> 188, RUE DU CADUCEE PARC EUROMEDECINE 34195 MONTPELLIER CEDEX 5 FRANCE</p> <p><b>PRIMO™</b> HOUSSE SUPERIEURE DETACHABLE 200 cm [79"] x 85 cm [33.5"]</p> <p><b>PRIMO™</b> DETACHABLE TOP COVER 200 cm [79"] x 85 cm [33.5"]</p> <p> = 250 Kg [551 lb]     = 250 Kg [551 lb]</p> <p> </p> <p></p> <p>E17170X0110P001</p>	<p><b>Hill-Rom Industries S.A</b> 188, RUE DU CADUCEE PARC EUROMEDECINE 34195 MONTPELLIER CEDEX 5 FRANCE</p> <p><b>PRIMO™</b> HOUSSE SUPERIEURE DETACHABLE 200 cm [79"] x 77 cm [30.5"]</p> <p><b>PRIMO™</b> DETACHABLE TOP COVER 200 cm [79"] x 77 cm [30.5"]</p> <p> = 250 Kg [551 lb]     = 250 Kg [551 lb]</p> <p> </p> <p></p> <p>E17185X0110P001</p>	<p><b>Hill-Rom Industries S.A</b> 188, RUE DU CADUCEE PARC EUROMEDECINE 34195 MONTPELLIER CEDEX 5 FRANCE</p> <p><b>PRIMO™</b> HOUSSE FEU SUPERIEURE DETACHABLE 200 cm [79"] x 85 cm [33.5"]</p> <p><b>PRIMO™</b> FIRE DETACHABLE COVER 200 cm [79"] x 85 cm [33.5"]</p> <p> = 250 Kg [551 lb]     = 250 Kg [551 lb]</p> <p> </p> <p></p> <p>E17184X0110P001</p>

Bottom Cover		
P02033/P02034	P02044	P02047
<p><b>Hill-Rom Industries S.A</b> 188, RUE DU CADUCEE PARC EUROMEDECINE 34195 MONTPELLIER CEDEX 5 FRANCE</p> <p><b>PRIMO™</b> HOUSSE INFERIEURE FERMEE 200 cm [79"] x 85 cm [33.5"]</p> <p><b>PRIMO™</b> CLOSED BOTTOM ENCASEMENT 200 cm [79"] x 85 cm [33.5"]</p> <p></p> <p></p> <p>E17215-1-0112P001</p>	<p><b>Hill-Rom Industries S.A</b> 188, RUE DU CADUCEE PARC EUROMEDECINE 34195 MONTPELLIER CEDEX 5 FRANCE</p> <p><b>PRIMO™</b> HOUSSE INFERIEURE FERMEE 200 cm [79"] x 77 cm [30.5"]</p> <p><b>PRIMO™</b> CLOSED BOTTOM ENCASEMENT 200 cm [79"] x 77 cm [30.5"]</p> <p></p> <p></p> <p>E17220X0110P001</p>	<p><b>Hill-Rom Industries S.A</b> 188, RUE DU CADUCEE PARC EUROMEDECINE 34195 MONTPELLIER CEDEX 5 FRANCE</p> <p><b>PRIMO™</b> HOUSSE INFERIEURE FEU FERMEE 200 cm [79"] x 85 cm [33.5"]</p> <p><b>PRIMO™</b> FIRE CLOSED BOTTOM ENCASEMENT 200 cm [79"] x 85 cm [33.5"]</p> <p></p> <p></p> <p>175186-1-0112P001</p>

See details of the symbols for cleaning and disinfection in section: "Disinfection, servicing" page 25.

# Installing the patient

## Before placing the patient on the Primo™

*Assess the various risks, including but not limited to the following (incomplete list containing risks related to reasonably foreseeable misuse):*

- *risk of entrapment;*
- *potential falls from the bed;*
- *patient in state of confusion;*
- *patient's learning ability;*
- *children (aged less than 12 or under 1.46 m tall);*
- *persons lacking the mental capacity to recognize risky actions;*
- *unauthorized persons.*

*Install and use for the first time in accordance with these instructions.*

## List of compatible bed frame

*The bed / mattress / siderails combination (and in particular their respective dimensions) must be examined in order to make sure that it meets the demands of the IEC 60601-2-52 standard and the "Hospital Bed Safety Workgroup" guide, and that the resulting combination does not alter the performance of the devices, their safety or their usability for the usage.*

## Accessories

The following accessories can be added to the Primo™ :

- Optional seat cushion;
- Quick inflate/deflate kit;
- Transport bag.

These accessories are described in chapter «Accessories and spare parts page 33».



## Installation

1. Unpack the control unit and the mattress.
2. Check that all the components are present and correct, and that the power cord is not damaged.
3. Place the rolled up mattress on the top of the bed at the foot end and unroll it.


**i** The CPR valve must be positioned at the foot end of the bed. The CPR valve must remain accessible for use in the event of emergency deflation. Check that the symbol  on the cover is at the foot of the bed.

**i** Keep the 2 elastic bands.

4. Attach the mattress using the straps at the head and foot ends.


5. Adjust the length of the straps to safely secure the mattress.



**Make sure that the mattress is correctly installed and securely attached, and in particular that it is well centered on the mattress support platform and pushed tightly against the foot board to avoid the danger of entrapment.**


**Make sure that the attachments of the mattress do not become entangled in the moving parts of the mattress support platform, such as actuators, CPR handles, etc. With articulated chassis, make sure that the straps of the mattress are only attached to the mobile head and foot sections and NOT to the main fixed chassis.**


6. Attach the control unit to the foot board of the bed (C) using the universal attachment support (D).

7. Open the cover to uncover the mattress.
8. Insert the hose integrated in the pocket (E) at the foot of the mattress.

9. Connect the air hose to the CPR valve (F). A click is heard when the hose is correctly connected.


10. Connect the power cord connector to the plug (G).
11. Turn the lock ring (H) clockwise to lock the coupling.
12. Check that the CPR valve is closed.

13. Make sure that the power cord is not damaged and connect it to the control unit.
14. Attach the power cord to immobile parts of the bed frame using the cable clamps. Make sure that the cord will not be jammed or drag along the floor.

***Check that the clamps of the power supply cord are attached to the main frame, and not to articulated parts, such as the head or foot sections. Failure to follow this recommendation could result in damage to the equipment.***

***Take care not to damage the power cord when moving the bed. Failure to follow this recommendation could result in corporal injury.***

15. Connect the power cord to the wall outlet.
16. Press the  button on the control unit to switch on the device.  
The yellow visual alarm signal  comes on.  
When connecting to the mains for the first time, the mattress automatically inflates. This takes about 15 minutes.

***Do not place the patient on the mattress while the indicator light is flashing during the initial inflation phase.***

***Placing a cotton sheet on the mattress will improve patient comfort and make it easier to administer care.***

***When the mattress is in use, make sure that the patient's sacrum is centered and positioned at the I-mmersion™ Therapy symbol located on each side of the mattress.***

# Mobilizing and securing the patient

---

## Functions

---

### Therapeutic mode

The Primo™ automatically works in continuous (or static) mode at low pressure.

The pressure in the air bladders is automatically regulated by the innovative I-mmersion™ pressure control system.

Irrespective of the positions of the bed frame articulations, this patented technology permanently detects the weight and the position of the patient and dynamically adjusts the pressure of the support accordingly.

The I-mmersion™ technology also includes an alarm that is activated whenever a fault occurs in the control system.

*The Primo™ is designed for optimal therapeutic benefits when the adjustable head section is inclined to between 0° and 45°.*

*The performance of the Primo™ diminishes at angles greater than 45°.*

**i** *Patients are ideally positioned when their hip is aligned with the patient position indicator on the bed.*

---

### Maximum inflation


This function increases the pressure inside the mattress in order to obtain a firm surface that facilitates patient movements when washing or changing position.

During the time of inflation, about 5 minutes, the effect of the I-mmersion™ sensor is canceled, and the inflation green indicator light flashes. The device automatically returns to therapy mode.

This function is also used to inflate the seat cushion (see “Seat cushion” page 33).


## CPR

In emergencies, opening the CPR valve (cardiopulmonary resuscitation) quickly deflates the mattress, providing a hard surface for external heart massages.

**Never allow a non-qualified person to operate this function and check that no obstacles (e.g., limbs, accessories, objects, power cords) or persons are under the head section.**

- i** The CPR valve is at the foot end of the mattress. It can be easily identified by its yellow plug and the symbol on the mattress cover.

1. If necessary, lower the siderails and the head section of the bed, or place the bed frame in the CPR position before using the emergency CPR function (refer to the instructions of the bed manufacturer).
2. Place a CPR board under the patient, or follow the protocol of the CPR function.
3. To activate the CPR, turn the yellow CPR valve clockwise as shown by the arrow (see opposite).

Activating the CPR quickly deflates the mattress.



- i** The yellow visual alarm signal  comes on to indicate the state of low pressure, and the audible alarm signal is activated.

### Cancel the CPR

1. Turn the CPR valve anticlockwise to its original position.

The device reverts to therapy mode.

2. If necessary, place the bed frame in the appropriate position (refer to the instructions for the bed frame).

## Alarms


### Alarm silence

 This control switches off the mains power fault and malfunction alarms.

Press this button to silence the alarms.

 Audible alarm is re-activated automatically after about 15 minutes.

The alarm can be stopped again for 15-minute periods until the problem is solved. The indicator light remains on. Refer to the Service Manual for more details.

Alarm silence

### Mains power fault


The visual alarm yellow  comes on if the unit is disconnected from the mains, in the event of a power failure, or during transport. A low-priority audible alarm sounds immediately, intermittently and at moderately strong volume.

 The button  remains active.

To stop the alarms, reconnect the control unit to the mains power supply, press  to silence the battery, or press  (see "Alarm silence" above).

### Malfunction

The yellow visual alarm  comes on after about 15 minutes if a pressure fault occurs. A low-priority audible alarm sounds immediately, intermittently and at moderately strong volume.

# Moving the device

## Transport mode

In order to move the patient on the bed, it is necessary to proceed with the following operations:

1. Press the inflate button to fully inflate the mattress. The indicator light near the  symbol flashes for 5 minutes time of inflation.
2. After the inflation finished, press  to switch off the control unit.
3. Disconnect the device from the mains by removing the plug from the wall socket.

**Never pull on the power cord or you may damage it. A damaged power cord is an electric shock hazard.**

**Make sure that the cord does not drag along the floor, do not run over it when moving the bed frame and take care not to trip up over it. If necessary, use the attachments provided with the device.**

The device automatically switches to transport mode and remains inflated, but the pressure is no longer regulated.

-  **The mattress is independant and can support the patient for 2 hours when disconnected from the mains.**

4. Safely stow the power cord.

**Make sure that the device does not remain disconnected from the mains for more than 2 hours when a patient is installed.**


5. Reconnect the device as soon as the destination is reached.
6. Press the  button. The device automatically returns to the normal operating mode.

-  **The CPR function remains available in transport mode. It is manual.**

## Transferring the mattress

Proceed as follows to transfer the device from one bed to another:

**Never leave the patient on the device.**

**Make sure that the bed frame brakes are applied to prevent any accidental movements.**

1. Press  to switch off the control unit, and to silence the battery.

-  **In this case no mains power fault alarm will be reported.**


2. Disconnect the power cord from the wall outlet.
3. Remove the cable clamps that attach the power cord to the bed.
4. Undo the two straps securing the mattress to the bed (at the head and foot ends).
5. Place the control unit on the mattress to prevent it from being damaged.
6. Transfer the mattress on the other bed.
7. To install the Primo™ on the other bed, refer to steps 4 to 16 of the installation procedure (page 16).

## Storing the mattress


***Make sure that the bed frame brakes are applied to prevent any accidental movements.***

1. Switch off the device by pressing  .  
The indicator light goes out.
2. Activate the CPR valve to help expel air from the mattress.
3. Disconnect the power cord from the wall outlet.

***Never pull on the power cord or you may damage it. A damaged power cord is an electric shock hazard.***

4. Remove the cable clamps that attach the power cord to the bed.
5. Disconnect the air hose and electric cable connector of the CPR valve.
6. Undo the two straps securing the mattress to the bed (at the head and foot ends).
7. Starting at the head end of the bed, roll up the mattress slowly to allow any remaining air to be expelled and place the rubber bands on each end of the rolled mattress.
8. Close the CPR valve to prepare the mattress for the next installation.
9. Store the mattress in its transport bag (see "Transport bag" page 36).

---

## Cleaning and Disinfection

---

### Safety Tips

- *Check that the brakes are applied on the bed on which the device is installed.*
- *Lock out all electrical functions.*
- *Disconnect the device and stow the power cord.*
- *Check that the connectors are securely connected to prevent water from entering the mattress.*
- *Never clean the device by pouring water on it, or with a high-pressure jet.*
- *Never use water at a temperature of more than 70°C.*
- *Avoid excess water on the connectors.*
- *Refer to the recommendations of the cleaning and disinfecting product manufacturer.*
- *Thoroughly dry the device before reusing it.*
- *Suitable individual protective equipment must be worn during the phases of the cleaning operations (blouse, gloves, eye protection, etc.).*

*Failure to implement one or more of these recommendations may lead to damage or deterioration, preventing use of the device and rendering the warranty void.*

---

### Infection control

**Insufficient cleaning = risk of infection (biological danger).**

All parts must be kept clean at all times in order to avoid the risk of infection. All necessary precautions must be taken to eliminate all visible soiling.

*The following instructions are not designed to replace more appropriate cleaning and disinfection protocols drawn up by the hygiene officer or by other bodies for your hospital in the event of particular infectious situations.*

## General recommendations

The cleaning and disinfecting method described below applies specifically to the device and its accessories and is designed to save time and to help combat nosocomial infection more effectively.

*Hill-Rom recommends that the Primo™ be disinfected before first use.*

*When cleaning, always check the cover for cuts, tears, cracks or snags. Never use a mattress with a damaged cover.*

*Using cleaning and disinfecting protocols or products other than those recommended by Hill-Rom can compromise the conformity of the devices and patient safety and render the warranty null and void.*

*Hill-Rom recommends that the Primo™ should be disinfected before scrapping, in accordance with applicable local regulations.*

**i** *Cleaning and disinfection are two separate procedures.*

### Products to be avoided:

Never use cleaning agents, detergents, degreasing agents or industrial solvents containing any of the following products in order to avoid damaging the device:

	Phenol		Hydrochloric, nitric or sulfuric acids		Dimethylformamide
	Cresol		Soda		Tetrahydrofuran

Do not use highly acidic detergents or disinfectants (pH<4).

Do not use highly basic detergents or disinfectants (pH>10).

Never use abrasive cleaning products or materials, such as scouring pads.

*The control unit is not protected against excessive ingress of fluids. Disconnect it before cleaning whenever using a water jet or there is a risk of damage.*


---

## Frequency of cleaning and disinfection of different parts of the mattress

Parts to be kept clean at all times.

+  - After a patient's departure.  
- If the inspection reveals ingress of fluids, every time the linen is changed.

---

## Clean and disinfect after the departure of a patient or when changing the linen.

### Recommended products for cleaning and/or disinfection:

**Terralin® Protect (Schülke)**

Concentration: 5 ml/l (= 0.5%), contact time: 1 h.

### Cleaning

Clean the device with a cloth slightly dampened with hot water and a neutral detergent solution. Check that the solution does not contain any of the products to be avoided mentioned above.  
(see "Products to be avoided:" page 26)

### Cleaning tough stains


**i** Quickly wipe away any traces of pharmaceutical solutions used for the patients to avoid damage to the surface.

- Remove tough stains using neutral detergents or a chlorine-based solution with a concentration less than or equal to 1,000 ppm and use a soft brush.
- To eliminate hardened stains (excreta, other forms of soiling), soften them by soaking and take care to thoroughly dry the cover before putting it back on the mattress.

## Disinfection

In the event of visible soiling, Hill-Rom recommends that the device should be disinfected with an intermediate disinfectant (tuberculocidal) that complies with the applicable regulations (e.g. with the requirements of directive 93/42/EEC).

For all other disinfectants, follow the recommendations below:

 <p>NaClO NaDCC C ≤ 1000 ppm</p>	Chlorine-based solutions can be used. Concentration must be less than or equal to 1,000 ppm.
 <p>C<sub>2</sub>H<sub>5</sub>OH 25% C<sub>2</sub>H<sub>5</sub>OH</p>	Ethanol-based solutions (alcohol) can be used. Concentration must not exceed ¼ of ethanol for ¾ of water.

**i** Proceed with the final steps before reusing the cleaned and disinfected mattress (see “Final steps” page 30).

---

## Cleaning and disinfection at regular intervals or in the event of high risk of contamination

Follow the same instructions as above (See “Clean and disinfect after the departure of a patient or when changing the linen.” page 27), but with the products listed below.

### Recommended products for cleaning and/or disinfection:

#### After the departure of patients with infectious diseases:

- Dismozon® Pur (**Bode**)  
Concentration: 40g/l/l (= 4%), contact time: 1 hour

#### After patients with Clostridium Difficile:

- Dismozon® Pur  
Concentration: 15g/l (= 1.5%), contact time: 2 hours
- **Anioxy-Spray WS (Anios)**  
Ready to use, contact time: 30 min

### Dry steam cleaning


The mattress can be dry steam cleaned.

**i** Dry steam or superheated steam contains no more than 6% of water in suspension and avoids the effects of condensation.

Nevertheless, to avoid damage due to high pressure or an abnormal surface temperature, take the following precautions:

- Use low pressure steam on the electric parts.

- Do not use accessories such as high pressure hoses (A). It is preferable to use soft non-metallic brushes (B) and wipers (C) in such a way as to reduce the pressure to an acceptable level.

- Prevent water and steam from getting into connectors that are not in use.
- Do not brush and apply reduced pressure to labels and markings.
- Dry carefully and look for signs of water ingress.
- Test the device before reusing it.

#### Hill-Rom recommends the use of the Sanivap® method:

- Application time recommended by the supplier: about 1 sec.
- Distance (recommended by the supplier) between the source and the surface to be cleaned: about 3 mm.
- Application: exterior parts of the mattress (covers).
- Frequency: for each new patient.

#### Machine washing the upper cover and the seat cushion

The upper cover of the mattress and the seat cushion can be machine washed. Follow the instructions below to prepare for machine washing.


*Do not always wash the upper cover and the cover of the seat cushion to avoid reducing its service life. The upper cover and seat cushion cover must only be machine washed if there is a particular risk of infection.*

	Maximum temperature: 70°C. Delicates.
	Only clean with oxygenated agents.
	Do not iron.
	Do not dry clean.
	Can be tumble dried. Moderate temperature: 50°C.


**i** Proceed with the final steps before reusing the cleaned and disinfected device.

---

## Final steps

*Always remove all traces of the products used when washing or disinfecting from the mattress.*


*Make sure that all the parts of the mattress are perfectly dry before installing to avoid any risk of condensation forming inside the mattress.*

---

## Maintenance

---

### Safety recommendations

*Never modify the device without Hill-Rom's prior written consent.*

*Only facility-authorized personnel should perform maintenance.*

*Changes made by unauthorized personnel may result in damage to the equipment and/or serious injury to staff or users.*

Before maintenance or servicing works:

- Check that the brakes are applied on the bed on which the device is installed.
- Lock out all electrical functions.
- Disconnect the device.
- Secure the mattress support platform and take whatever steps are necessary to prevent any movement.

For the complete maintenance procedure, refer to the Primo™ Service Manual. Contact your local Hill-Rom representative (see back of manual).

---

### Preventive maintenance

*A Service Manual and spare parts catalog are provided with the device. The wiring diagrams, list of components, descriptions and calibration instructions can be obtained on demand from Hill-Rom After-Sales.*

*The frequency of inspections must be adapted to the general condition of the product and its use, for example, if the device is used by heavy patients. It is the responsibility of the facility to implement a preventive maintenance program for the device according to its conditions of use.*

The device and accessories should be inspected at least once a year to keep them in good condition and working properly.

Every three years, it is preferable to ask Hill-Rom After-Sales Service or a Hill-Rom approved supplier to inspect the device in order to keep it in safe and good working order over time. Depending on the maintenance operations and observations, the date of the next inspection must be recommended every time the bed is serviced.


In order to benefit from optimal and rapid service when calling Hill-Rom about your Primo™, provide the serial number of the Primo™ about which you are calling (See “Product identification labels” page 13).

In view of the specific conditions of use in domestic care environments, the expected life time of the device, its parts and the accessories is 3 years.

---

## Troubleshooting

The Primo™ is designed to function automatically. Consequently, device troubleshooting is limited to a few checks.


***Always disconnect the device before troubleshooting.***

In the event of a malfunction, check the following points:

- the electrical connection with the control unit;
- the power cord is not damaged;
- the pneumatic and electric connections between the mattress and the control unit are OK;
- the CPR valve is properly closed;
- the mattress is not damaged (tears or cuts);
- the air filter is clean. Clean the air filter every 6 months.
- the fuse is OK.

**i** Refer to the Primo™ Service Manual for detailed troubleshooting instructions.

### Cleaning and replacing the air filter


The reusable air filter can be washed with water and a cleaning agent.

Thoroughly wash and dry the air filter before putting it back in place.

**i** The air filter can be ordered separately. For more information, refer to the list of spare parts of the Primo™.

## Replacing the fuse on the power socket

1. Disconnect the device.
  2. Remove the power cord of the connector.
  3. Remove the fuse box located on the power socket using a small screwdriver.
- i** Check that the new fuse matches the characteristics on the label and meets the IEC 60269-1 standard.
4. Replace the old fuse.
  5. Securely close the fuse box.
- i** Two others fuses are located on the alimentation card. Refer to the Primo™ Service Manual for replacement.

If the problem is not solved, contact your local Hill-Rom representative (see back of manual). Indicate the serial number of the device, which can be found on the control unit, when requesting repairs.

---

## De-commissioning

The device and its accessories should be cleaned and disinfected before de-commissioning.

Decommissioned equipment materials (plastics, electrical components, etc.) must be recycled in accordance with local recycling regulations. Always meet the applicable demands and local rules relating to environmental protection, in particular for waste from medical devices.

Do not dispose of electric and electronic equipment in the waste bin (as per directive 2012/19/EU).

Never discard the batteries or accumulators of your device. They may contain substances and metals that are hazardous for the environment and health (as per Directive 2006/96/EC).

- i** Refer to the Primo™ Service Manual for detailed battery replacement instructions.


# Accessories and spare parts

## Seat cushion

The optional seat cushion is an accessory that can be inflated using the control unit of the Primo™ mattress.


### Installation

1. Place the seat cushion at the foot of the bed.
2. Open the Velcro® straps (A) on the seat cushion, disconnect and stretch the spiral hose (B).
3. Connect the cushion's spiral hose to the control unit. The hose clicks into place. 
4. Connect the control unit to the mains.
5. Press the button to inflate the cushion on the control unit. The cushion inflates. The inflation indicator light automatically flashes for 5 minutes, then stops. 
6. Once the cushion is firmly inflated, press the cushion connector lock (C) and disconnect the coupling.
7. Place the seat cushion on the patient's chair and secure using the strap provided (D).

Patients can now sit on the cushion, which makes a whistling sound. The cushion is optimally inflated when the whistling stops.

8. Connect the cushion hose to the non-return valve (E) of the seat cushion.


9. Wind the Velcro® straps (A) round the tube and close.

---

## Deflating and packing the seat cushion

1. Open the cushion cover by separating the Velcro® straps on the folds on the side of the cushion.
2. Locate the blue deflation valve and loosen the plug (G) by turning through a ½ turn anticlockwise.

3. Starting on the side opposite the valve, fold the cushion in order to expel the air.
4. Close the plug and attach the Velcro® straps.
5. Store the seat cushion.
6. The device must be stored in its original packaging:
  - protected against light and damp;
  - at least 10 cm above floor level to prevent fluid ingress;
  - protected against dust;
  - outside passageways.

---

## Cleaning and Disinfection


*When cleaning, always check the seat cushion cover for cuts, tears, cracks or snags. Hill-Rom recommends that damaged underlays should be scrapped.*

Refer to chapter "Disinfection, servicing" page 25, up to paragraph "Disinfection" page 28.


## Upper covers

Upper covers for each mattress model

Mattress model	Upper cover
P02033B	Upper cover, 200*85*16 cm
P02034B	Upper cover, 200*85*16 cm
P02044B	Upper cover, 200*77*16 cm
P02047B	Upper cover, fire, 200*85*16 cm

## Quick inflate/deflate kit


Refer to the 170196 User Manual for the use of the quick inflate/deflate kit.

---

## Transport bag

You can order a transport and storage bag for the Primo™.

- i** *All the accessories can be ordered separately.  
Refer to the list of spares for the product part numbers.  
Contact your local Hill-Rom representative or your local reseller (see back of manual).*

---

## Warranty and after sales service conditions

The warranty for our products will be rendered null and void, in part or in total, in the event of:

- repairs, installations, assembly, modifications or checks and tests are not conducted by the manufacturer's maintenance personnel or by personnel authorized by the manufacturer;
- the electric system does not provide the conditions allowing for the use of devices that comply with the standard EN 60601-1. In hospital wards in particular, the electric system does not meet the requirements applying to medical establishments;
- the device is not used in accordance with this manual;
- accessories are used that do not meet the requirements of this manual.

The details of after sales service contacts in all countries are shown on the back of this manual.

---

## Regulatory requirements

**CE** The Primo™ is a class IIa medical device that complies with the requirements of Directive 93/42/EEC and its amendments.

0459 The CE mark was applied in 1998.

The Primo™ is designed and manufactured in accordance with the following standards and classifications:

Quality standards:

- ISO 9001 : 2008                      - ISO 13485 : 2012                      - ISO 14001 : 2004

Technical standards:

<b>Name</b>	<b>P02033/P02034/P02044</b>	<b>P02047</b>
EN IEC 60601-1: 2006	YES	YES
EN IEC 60601-1-2: 2007	YES	YES
EN IEC 60601-1-6: 2010	YES	YES
EN IEC 60601-1-8: 2007	YES	YES
EN IEC 60601-1-11: 2010	YES	YES
EN ISO 14971: 2012	YES	YES
EN ISO 10993-1: 2009	YES	YES
EN ISO 10993-5: 2009	YES	YES
EN ISO 10993-10: 2009	YES	YES
EN ISO 15223-1: 2012	YES	YES
EN 597-1: 1995	YES	YES
EN 597-2: 1995	YES	YES
SS 876 00 01 (NT FIRE 037)	YES	YES
UNI 9175: 2008 (Class 1 IM)	NO	YES
BS 6807: 2006, Clause 9, Ignition/5	NO	YES
BS 7177: 1996 (Medium hazard)	NO	YES

# Manufacturer's declaration - Electromagnetic Compatibility

## Electromagnetic emissions compliance

Manufacturer's recommendations and declaration – electromagnetic emissions		
The Primo™ is intended for use in the electromagnetic environment specified below. The customer or the user of the Primo™ should make sure it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment - Recommendations
RF emissions CISPR 11	Group 1	The Primo™ only uses RF energy internally. Consequently, it only produces very weak RF emissions that are unlikely to cause interference with nearby electronic equipment.
RF emissions CISPR 11	Class B	The Primo™ is suitable for use in all establishments, including domestic premises and premises directly connected to the public low-voltage power supply in residential buildings.
Harmonic emissions IEC 61000-3-2	Class A	
Emissions of voltage fluctuations/flicker IEC 61000-3-3	Compliant	

## Electromagnetic immunity compliance

Manufacturer's recommendations and declaration – electromagnetic immunity			
The Primo™ is intended for use in the electromagnetic environment specified below. The customer or the user of the Primo™ should make sure it is used in such an environment.			
Immunity test	Test level IEC 60601	Compliance	Electromagnetic environment - Recommendations
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV on contact ±8 kV in the air	±6 kV on contact ±8 kV in the air	Floors should be made of wood, concrete or ceramic tiles. If the floor is covered with a synthetic material, relative humidity must be at least 30%.
Fast electric transients in bursts IEC 61000-4-4	±2 kV on power supply lines ±1 kV for the input/output lines	±2 kV on power supply lines ±1 kV for the input/output lines	The quality of the electric power supply should be the equivalent of that found in commercial or hospital environments.
Shock waves IEC 61000-4-5	±1 kV between line(s) ±2 kV between line(s) and earth	±1 kV between line(s) ±2 kV between line(s) and earth	The quality of the electric power supply should be the equivalent of that found in commercial or hospital environments.
Magnetic field at the frequency of the mains power supply (50/60 Hz) IEC 61000-4-8	3 A/m	3 A/m	The magnetic field at the frequency of the mains supply must be of a level that is characteristic of a place located in a commercial or hospital environment.

Manufacturer's recommendations and declaration – electromagnetic immunity			
The Primo™ is intended for use in the electromagnetic environment specified below. The customer or the user of the Primo™ should make sure it is used in such an environment.			
Immunity test	Test level IEC 60601	Compliance	Electromagnetic environment - Recommendations
Voltage dips, brief outages and variations in voltage on the electric power supply lines IEC 61000-4-11	<ul style="list-style-type: none"> <li>• &lt;5% <math>U_T^a</math> (dips &gt;95% of <math>U_T</math>) for 0.5 cycle</li> <li>• 40% <math>U_T</math> (dips = 60% of <math>U_T</math>) for 5 cycles</li> <li>• 70% <math>U_T</math> (dips = 30% of <math>U_T</math>) for 25 cycles</li> <li>• &lt;5% <math>U_T</math> (dips &gt;95% of <math>U_T</math>) for 5 seconds</li> </ul>	<ul style="list-style-type: none"> <li>• &lt;5% <math>U_T</math> (dips &gt;95% of <math>U_T</math>) for 0.5 cycle</li> <li>• 40% <math>U_T</math> (dips = 60% of <math>U_T</math>) for 5 cycles</li> <li>• 70% <math>U_T</math> (dips = 30% of <math>U_T</math>) for 25 cycles</li> <li>• &lt;5% <math>U_T</math> (dips &gt;95% of <math>U_T</math>) for 5 seconds</li> </ul>	The quality of the electric power supply should be the equivalent of that found in commercial or hospital environments. If the user of the Primo™ requires the appliance to be in continuous use during power cuts, it is advisable to power the Primo™ from an uninterruptible power source or a battery.

a.  $U_T$  is the nominal value of the supply voltage applied during the test.


## Recommended separation distances

RF communication equipment (including peripherals such as antenna cables and external antennas) should be used no closer to any part of the Primo™, including cables specified by Hill-Rom, than the recommended separation distances shown in the following table. Otherwise, degradation of the performance of this equipment could result.

Wireless communication system	Separation distance (m)
TETRA 400 (limited to 10 W ERP*)	0,3
Public safety(460-470 MHz), GMRS 460 (limited to 5 W ERP*)	0,2
GMRS 460 (limitede to 2 W ERP*)	0,1
GSM 850, GSM 900, RFID 868 MHz, TETRA 800 (limited to 2 W ERP*)	0,3
GSM 1900 (limited to 1 W ERP*)	0,2
WLAN 802.11a 5 GHz (limited to 1W ERP*)	0,7
iDEN 820, CDMA 850, GSM 1800, CDMA 1900 (limited to 0,6 W ERP*)	0,2
FRS 460 (limited to 0,6 W ERP*), PMR 446	0,1
UMTS, DECT (limited to 0.25 W ERP*)	0,1

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© 2012 Hill-Rom Services, Inc. TOUS DROITS RESERVES.

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