



# T Frame Gantry

## Technical Specification

Introduced in 1995 the long serving and reliable T Frame Gantry has stood the test of time with its ease of installation and maintenance it is always a popular choice.

The T Frame Gantry combines with the Wispa Lite to create a temporary or emergency solution, or where a fixed ceiling track installation or mobile hoist is not suitable. A lightweight and portable moving and handling solution for any care environment.

The T Frame Gantry has a safe working load of 200kg (31 stone) and can be adjusted in height between 2060mm and 2380mm

### Easy Assembly Easy Use

Can be assembled or disassembled by one person within five minutes

Operates with a battery powered hoist.

### Versatile


Suits the majority of clients  
The T Frame gantry system complements the Wispa Lite portable hoist and patient slings with its maximum load of 200kg (31 stones) (440lb).

### Installation

No installation requirement, electrical connection, or fixings and can be removed in minutes.

Totally self supporting and independent of the building structure.

Provides an alternative where  
Ceiling construction limits fixed track installation and user preference.



### Temporary Solution

Ideal temporary solution  
No power supply

### Affordable

An Affordable alternative to permanent track installation


### Support Services

Free of charge loan store \ OT training and demonstration visits are available upon request.

Service and maintenance programmes are available.

# T Frame Gantry

## Technical Specification



A	Overall Width	2,580mm
B	Internal Width (max)	2,500mm
C	Lowest Height	2,060mm
D	Highest Height	2,380mm
E	Depth	1,120mm